



## Commonly Used Frontier Consumables



### ITF Needle Set N

Deactivated interface needle for PY-3030, PY-2020, and Rx-3050. Set of 3.

**SKU:** PY1-1274



### Glass Insert

Surface deactivated, with deactivated quartz wool.


**SKU:** PY1-3337



### Quartz Pyrolysis Tube

For pyrolyzer models PY-3030D and PY-3030S, 120 mm (L), made of quartz.


**SKU:** PY1-3018A



### O-Ring P-6R

Used to seal quartz pyrolysis tubes. Sold in a pack of 10.

**SKU:** PY1-2021



### Graphite Vespel Ferrule

6 mm hole i.d., for interface union. Includes 5 pieces

**SKU:** PY1-7911



### Ultra ALLOY Capillary Column

30 meters, 0.25 mm i.d., 0.25  $\mu$ m, 1 each. 5% diphenyl 95% dimethylpolysiloxane

**SKU:** UA5-30M-0.25F



### Eco Cup LF

Deactivated stainless steel sample cup (80  $\mu$ L). Bottle contains 100 pieces.

**SKU:** PY1-EC80F



### Eco-Stick DF

For stainless steel Eco-Cup LF. Includes 50 pieces.


**SKU:** PY1-ES20F



### Vent-free GC/MS Adapter

Vent-Free GC/MS Adapter N-50, recommended for Agilent/Thermo GC.

**SKU:** MS402280



### VF Metal Ferrule D

Used to connect the vent-free tube with a separation column (or an EGA tube). 20 pieces.

**SKU:** MS402167



### Ultra ALLOY EGA Capillary Tube

UA deactivated metal capillary tube, 0.15 mm i.d. x 0.47 mm o.d., 2.5 meters.

**SKU:** UADTM-2.5N



### SS Adapter Packing

Selective Sampler adapter packing (seal). Comes with 10 each.

**SKU:** SS1-1011



**Quantum**  
Analytics

Visit online store for full Frontier consumables catalog

**store.lqa.com**